1.09204.0100
1.09204.0103
1.09204.0500
1.09204.0503
1.09204.1000
1.09204.1022
1.09204.2500
1.09204.9025

Microscopy

Giemsa's azur eosin methylene blue solution

for microscopy

For professional use only

IVD *In vitro* Diagnostic Medical Device

Intended purpose

This "Giemsa's azur eosin methylene blue solution - for microscopy" is used for human-medical cell diagnosis and serves the hematological, clinico-cytological, and histological investigation of sample material of human origin. It is a staining solution that when used together with other *in vitro* diagnostic products from our portfolio makes target structures evaluable for diagnostic purposes (by fixing, embedding, staining, counterstaining, mounting) in human-hematological, histological and clinico-cytological specimen materials, for example smears of whole blood and bone marrow, as well as paraffin sections.

Unstained structures are relatively low in contrast and are extremely difficult to distinguish under the light microscope. The images created using the staining solutions help the authorized and qualified investigator to better define the form and structure in such cases. Further tests must be carried out according to recognized, valid methods to reach a definitive diagnosis.

Principle

When used in hematological applications, Giemsa's stain is frequently used in combination with other staining solutions, e.g. with May-Grünwald's solution for Pappenheim (MGG) overview staining. This staining solution generally stains the nuclei red, based on the molecular interaction between the eosin Y dye and an Azure B-DNA complex. Both dyes assemble to an Eosin Y - Azure B-DNA complex and the intensity of the resulting stain depends on the content of Azure B and the ratio of Azure B : Eosin Y. Furthermore, the resulting stain can vary depending on the influence of fixation, staining times, pH-value of the solutions or buffer substances. In histology and clinico-cytological applications, Giemsa's staining without additional stainings is used as an extended overview staining method. In this method, the color of the various cell components is influenced by pretreatment of the specimen material. Here, chromatin-containing structures (e.g. cell nuclei) appear in various blue shades, while the acidophilic components show up in a variety of red shades.

Sample material

Sections of formalin fixed, paraffin embedded tissue (3 - 4 µm thick paraffin sections) or fresh, native whole blood and bone-marrow smears as well as clinical cytological material like urine sediment, sputum, smears from fine needle aspiration biopsies (FNAB), rinses, imprints are used as starting material.

Reagents

Cat. No. 1.09204
Giemsa's azur eosin methylene blue solution 100 ml, 500 ml, 1 l, 2.5 l
for microscopy

Also required:

Cat. No. 1.06009 Methanol for analysis EMSURE® 1 l, 2.5 l, 5 l
ACS,ISO,Reag. Ph Eur
Cat. No. 1.09468 Buffer tablets pH 7.2 100 tabs
for preparing buffer solution acc. to WEISE
for staining of blood smears

or

Cat. No. 1.11373 Buffer tablets pH 6.4 100 tabs
for preparing buffer solution acc. to WEISE
for the staining of blood smears

or

Cat. No. 1.11374 Buffer tablets pH 6.8 100 tabs
for preparing buffer solution acc. to WEISE for
the staining of blood smears

for the staining of paraffin sections:

Cat. No. 1.00063 Acetic acid (glacial) 100% anhydrous 1 l, 2.5 l
for analysis EMSURE®
ACS,ISO,Reag. Ph Eur

Cat. No. 1.09634 2-Propanol for analysis EMSURE® 1 l, 2.5 l, 5 l
ACS,ISO,Reag. Ph Eur

for the Pappenheim staining:

Cat. No. 1.01424 May-Grünwald's eosine-methylene 100 ml, 500 ml,
blue solution modified 1 l, 2.5 l
for microscopy

Sample preparation

The sampling must be performed by qualified personnel.

All samples must be treated using state-of-the-art technology.

All samples must be clearly labeled.

Suitable instruments must be used for taking samples and their preparation. Follow the manufacturer's instructions for application / use.

When using the corresponding auxiliary reagents, the corresponding instructions for use must be observed.

Deparaffinize and rehydrate paraffin sections in the conventional manner.

Notes on Giemsa staining of paraffin sections

Always employ separate xylene or Neo-Clear™ (Cat. No. 1.09843) rinse baths when Giemsa staining paraffin sections as any ethanol traces in the solutions may result in the preparations being discolored.

Pretreatment of bone marrow and iliac crest biopsy materials

Optimal results can be achieved using an OSTEOSOFT® mild decalcifier-solution (Cat. No. 1.01728).

To gently remove any calcification, the fixed biopsy materials are first placed in OSTEOSOFT® for 18 - 24 hours, after which they are transferred to histoprocessing. Blocks are carefully cut and, if required, are again treated with OSTEOSOFT® for an additional 20 minutes.

Reagent preparation

Giemsa's azur eosin methylene blue solution

The solution is supplied as a concentrated staining solution and before use must be diluted with a buffer solution as described below. The diluted staining solution should be filtered before use.

Buffer solution

For preparation of approx. 1000 ml of solution, add and dissolve:

Buffer tablet, Cat. No. 1.11373 (pH 6.4), Cat. No. 1.11374 (pH 6.8) or Cat. No. 1.09468 (pH 7.2) depending on the required reaction color	1 tablet
Distilled water	1000 ml

Dilute Giemsa's staining solution for manual staining

For preparation of approx. 200 ml solution mix:

Giemsa's azur eosin methylene blue solution	10 ml
Buffer solution	190 ml
Mix well, leave to stand for 10 min, and filter if necessary	

Dilute Giemsa's staining solution for staining with automatic stainer

For preparation of approx. 300 ml solution mix:

Giemsa's azur eosin methylene blue solution	25 ml
Buffer solution	275 ml
Mix well, leave to stand for 10 min, and filter if necessary	

In many cases precipitates of the dye form in the diluted staining solutions; these can be eliminated by repeating the filtration process.

Acetic acid 0.1 %, aqueous

For preparation of approx. 1000 ml solution mix:

Acetic acid 100 %	1 ml
Distilled water	1000 ml

Giemsa staining

Procedure

Air-dried smears

Staining in the staining cell / on the staining rack

The slides should be allowed to drip off well after the individual staining steps, as a measure to avoid any unnecessary cross-contamination of solutions.

The stated times should be adhered to in order to guarantee an optimal staining result.

Slide with air-dried smear	
Methanol	3 min
Dilute Giemsa's staining solution for manual staining	20 min
Buffer solution	1 min
Buffer solution	1 min
Air-dry (e.g. over night or at 50°C in the drying cabinet)	

Staining in the automatic stainer

The stated times should be adhered to in order to guarantee an optimal staining result.

	Time	Station	Dip
Slide with air-dried smear			
Methanol	3 min	2	on
Dilute Giemsa's staining solution for staining with automatic stainer	20 min	3	on
Buffer solution	1 min	4	on
Running tap water	2 min	5	on
Dry	3 min	6	-

All diluted solutions should be replaced after each working day.

Covering with non-aqueous mounting media (e.g. Neo-Mount™, DPX new, or Entellan™ new) and a cover glass is recommended for the storage of hematological specimens for several months.

After dehydration (ascending alcohol series) and clarification with xylene or Neo-Clear™, cytological samples can be mounted with water-free mounting agents (e.g. Entellan™ new, DPX new, or Neo-Mount™) and a cover glass and can then be stored.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

Result

	Buffer solution pH 6.4	Buffer solution pH 6.8	Buffer solution pH 7.2
Nuclei resp. chromatin	red to violet	red to violet	red to violet
Cytoplasm of lymphocytes	blue	blue	blue
Cytoplasm of monocytes	grey-blue	grey-blue	grey-blue
Neutrophilic granules	light violet	light violet	light violet
Eosinophilic granules	reddish to red-brown	reddish to red-brown	reddish to red-brown
Basophilic granules	dark violet	dark violet	dark violet
Thrombocytes	violet	violet	violet
Erythrocytes	reddish	reddish	reddish-brownish

Pappenheim staining

with May-Grünwald's solution and Giemsa's solution

Procedure

Air-dried smears

Staining in the staining cell

The slides must be immersed and moved briefly in the solutions, simple immersion alone yields inadequate staining results.

The slides should be allowed to drip off well after the individual staining steps, as a measure to avoid any unnecessary cross-contamination of solutions.

The stated times should be adhered to in order to guarantee an optimal staining result.

Slide with air-dried smear	
May-Grünwald's eosine-methylene blue solution modified	3 min
Dilute Giemsa's staining solution for manual staining	20 min
Buffer solution	1 min
Buffer solution	1 min
Air-dry (e.g. over night or at 50°C in the drying cabinet)	

Staining on the staining rack

The stated times should be adhered to in order to guarantee an optimal staining result.

Slide with air-dried smear			
May-Grünwald's eosine-methylene blue solution modified		cover completely	3 min
Buffer solution	1 ml	mix	
Dilute Giemsa's staining solution for manual staining		cover completely	20 min
Buffer solution		rinse	
Air-dry (e.g. over night or at 50°C in the drying cabinet)			

All diluted solutions should be replaced after each working day. Only in the case of May-Grünwald's eosin methylene blue solution modified should the concentrated solution, when used daily, be renewed after one working week, at the latest, or otherwise when necessary. Concentrated May-Grünwald's eosin methylene blue solution modified (in the event of evaporation) may not be refilled, otherwise the concentration of the staining solution is no longer correct.

Covering with non-aqueous mounting media (e.g. Neo-Mount™, DPX new, or Entellan™ new) and a cover glass is recommended for the storage of hematological specimens for several months.

After dehydration (ascending alcohol series) and clarification with xylene or Neo-Clear™, cytological samples can be mounted with water-free mounting agents (e.g. Entellan™ new, DPX new, or Neo-Mount™) and a cover glass and can then be stored.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

Result

	Buffer solution pH 6.4	Buffer solution pH 6.8	Buffer solution pH 7.2
Nuclei resp. chromatin	red-violet	purple to violet	violet
Cytoplasm of lymphocytes	blue	blue	blue
Cytoplasm of monocytes	grey-blue	grey-blue	grey-blue
Neutrophilic granules	light violet	light violet	violet
Eosinophilic granules	brick-red	brick-red	red-brown
Basophilic granules	dark violet	dark violet to black	dark violet to black
Thrombocytes	violet	violet	violet
Erythrocytes	reddish	reddish	reddish-grey

Giemsa staining

Procedure

Paraffin sections of punched iliac crest specimens and detection of *Helicobacter pylori*

Staining in the staining cell

The slides should be allowed to drip off well after the individual staining steps, as a measure to avoid any unnecessary cross-contamination of solutions.

Deparaffinize histological slides in the conventional manner and rehydrate in a descending alcohol series.

The stated times should be adhered to in order to guarantee an optimal staining result.

Always employ separate xylene or Neo-Clear™ (Cat. No. 1.09843) rinse baths when Giemsa staining paraffin sections as any ethanol traces in the solutions may result in the preparations being discolored.

Slide with histological specimen	
Distilled water	10 sec
Giemsa's azur eosin methylene blue solution (undiluted, filtered)	15 min
Acetic acid 0.1 %	10 sec
Distilled water	10 sec
2-Propanol	10 sec
2-Propanol	10 sec
2-Propanol	10 sec
Xylene or Neo-Clear™	5 min
Xylene or Neo-Clear™	5 min
Mount the Neo-Clear™-wet slides with Neo-Mount™ or the xylene-wet slides with e. g. Entellan™ new and cover glass.	

In the case of Giemsa's azur eosin methylene blue solution should the concentrated solution, when used daily, be renewed after one working week, at the latest, or otherwise when necessary. Concentrated Giemsa's azur eosin methylene blue solution (in the event of evaporation) may not be refilled, otherwise the concentration of the staining solution is no longer correct.

After dehydration (ascending alcohol series) and clarification with xylene or Neo-Clear™, histological slides can be covered with non-aqueous mounting agents (e. g. Entellan™ new, Neo-Mount™) and a cover glass and can then be stored.

The use of immersion oil is recommended for the analysis of stained slides with a microscopic magnification >40x.

Result

Nuclei, cells	blue to dark blue
Collagen, osteoid	pale blue
Eosinophilic granules	red
Acidic mucopolysaccharides, mastocytes granules, cartilage matrix	reddish-violet
Acidophilic materials	orange-red
<i>Helicobacter pylori</i>	blue to dark blue

Technical notes

The microscope used should meet the requirements of a medical diagnostic laboratory.

When using histoprocessors and automatic staining systems, please follow the instructions for use supplied by the supplier of the system and software. The diluted staining solution should be filtered before use.

Remove surplus immersion oil before filing.

Analytical performance characteristics

"Giemsa's azur eosin methylene blue solution" stains and thereby visualizes biological structures, as described in the "Result" chapters of this IFU. The use of the product is only to be carried out by authorized and qualified persons, this includes, among other things, sample and reagent preparation, sample handling, histoprocessing, decisions regarding suitable controls and more.

The analytical performance of the product is confirmed by testing each production batch. The successful participation in international interlaboratory tests on a regular basis provide an additional and unaffiliated confirmation of analytical specificity and repeatability.

For the following stains, the analytical performance was confirmed in terms of specificity, sensitivity and repeatability of the product with a rate of 100 %:

	Inter-assay Specificity	Inter-assay Sensitivity	Intra-assay Specificity	Intra-assay Sensitivity
Hematological staining				
Erythrocytes	20/20	20/20	14/14	14/14
Nuclei	20/20	20/20	14/14	14/14
Eosinophilic granules	20/20	20/20	14/14	14/14
Neutrophilic granules	20/20	20/20	14/14	14/14
Lymphocyte cytoplasm	20/20	20/20	14/14	14/14
Monocyte cytoplasm	20/20	20/20	14/14	14/14
Thrombocytes	20/20	20/20	14/14	14/14
Basophilic granules	20/20	20/20	14/14	14/14
Blood parasite nuclei	20/20	20/20	14/14	14/14
Histological staining				
Nuclei	20/20	20/20	14/14	14/14
Cells	20/20	20/20	14/14	14/14
Collagen	20/20	20/20	14/14	14/14
Osteoid	20/20	20/20	14/14	14/14
Eosinophilic granules	20/20	20/20	14/14	14/14
Acidic mucopolysaccharides	20/20	20/20	14/14	14/14
Mast cell granules	20/20	20/20	14/14	14/14
Cartilage matrix	20/20	20/20	14/14	14/14
Acidophilic materials	20/20	20/20	14/14	14/14
<i>Helicobacter pylori</i>	20/20	20/20	14/14	14/14

Analytical performance results

Intra- (performed on the same batch) and inter-assay (performed on different batches) data list the number of correctly stained structures in relation to the number of performed assays.

Clinical performance characteristics

Furthermore, the clinical performance of this product has been successfully proven in multiple scientific publications.

The diagnostic interpretation of the staining results, however, is to be carried out by qualified and authorized professionals, taking into account patient anamnesis, morphology, the use of adequate controls, and additional diagnostic tests, if appropriate. This method can be supplementarily used in human diagnostics.

The clinical performance of the use of the Giemsa's azur eosin methylene blue solution for the detection of *Helicobacter pylori*, for example, was compared with the immunohistochemical detection of *H. pylori*, since no "gold standard" for its diagnosis exists. The sensitivity and specificity of the Giemsa's azur eosin methylene blue solution for the detection of *H. pylori* in comparison to the immunohistochemical detection was determined as follows:

	Giemsa staining	Immunohistochemistry
Sensitivity	13/15	15/15
Specificity	15/15	15/15

Sensitivity: 13 samples out of 15: 86.7 %

Specificity: 15 samples out of 15: 100 %

Positive predictive value (PPV): 100 %

Negative predictive value (NPV): 88.3 %

The results of this Performance Evaluation confirms that the product is suitable for the intended use and performs reliably.

Diagnostics

Diagnoses are to be made only by authorized and qualified personnel. Valid nomenclatures must be used. This method can be supplementarily used in human diagnostics. Further tests must be selected and implemented according to recognized methods. Suitable controls should be conducted with each application in order to avoid an incorrect result.

Storage

Store the Giemsa's azur eosin methylene blue solution - for microscopy at +15 °C to +25 °C.

Shelf-life

The Giemsa's azur eosin methylene blue solution - for microscopy can be used until the stated expiry date.

After first opening of the bottle, the contents can be used up to the stated expiry date when stored at +15 °C to +25 °C.

The bottles must be kept tightly closed at all times.

Capacity

3500 - 5000 stainings / 500 ml

Additional instructions

For professional use only.

In order to avoid errors, the application must be carried out by qualified personnel only. National guidelines for work safety and quality assurance must be followed. Microscopes equipped according to the standard must be used.

Protection against infection

Effective measures must be taken to protect against infection in line with laboratory guidelines.

Instructions for disposal

The package must be disposed of in accordance with the current disposal guidelines.

Used solutions and solutions that are past their shelf-life must be disposed of as special waste in accordance with local guidelines. Information on disposal can be obtained under the Quick Link "Hints for Disposal of Microscopy Products" at www.microscopy-products.com. Within the EU the currently applicable REGULATION (EC) No 1272/2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 applies.

Auxiliary reagents

Cat. No. 1.00063	Acetic acid (glacial) 100% anhydrous for analysis EMSURE® ACS,ISO,Reag. Ph Eur	1 l, 2.5 l
Cat. No. 1.00496	Formaldehyde solution 4%, buffered, pH 6.9 (approx. 10% Formalin solution) for histology	350 ml and 700 ml (in bottle with wide neck), 5 l, 10 l, 10 l Titripac®
Cat. No. 1.00579	DPX new non-aqueous mounting medium for microscopy	500 ml
Cat. No. 1.00974	Ethanol denatured with about 1 % methyl ethyl ketone for analysis EMSURE®	1 l, 2.5 l
Cat. No. 1.01424	May-Grünwald's eosine-methylene blue solution modified for microscopy	100 ml, 500 ml, 1 l, 2.5 l
Cat. No. 1.01728	OSTEOSOFT® mild decalcifier-solution for histology	1 l, 10 l Titripac®
Cat. No. 1.03699	Immersion oil Type N acc. to ISO 8036 for microscopy	100-ml dropping bottle
Cat. No. 1.04699	Immersion oil for microscopy	100-ml dropping bottle, 100 ml, 500 ml
Cat. No. 1.06009	Methanol for analysis EMSURE® ACS,ISO,Reag. Ph Eur	1 l, 2.5 l, 5 l
Cat. No. 1.07961	Entellan™ new rapid mounting medium for microscopy	100 ml, 500 ml, 1 l
Cat. No. 1.08298	Xylene (isomeric mixture) for histology	4 l
Cat. No. 1.09016	Neo-Mount™ anhydrous mounting medium for microscopy	100-ml dropping bottle, 500 ml

Cat. No. 1.09203	Giemsa's azur-eosin-methylene blue for microscopy	25 g, 100 g
Cat. No. 1.09468	Buffer tablets pH 7.2 for preparing buffer solution acc. to WEISE for staining of blood smears	100 tabs
Cat. No. 1.09634	2-Propanol for analysis EMSURE® ACS, ISO, Reag. Ph Eur	1 l, 2.5 l, 5 l
Cat. No. 1.09843	Neo-Clear™ (xylene substitute) for microscopy	5 l
Cat. No. 1.11373	Buffer tablets pH 6.4 for preparing buffer solution acc. to WEISE for the staining of blood smears	100 tabs
Cat. No. 1.11374	Buffer tablets pH 6.8 for preparing buffer solution acc. to WEISE for the staining of blood smears	100 tabs
Cat. No. 1.11609	Histosec™ pastilles solidification point 56-58°C embedding agent for histology	1 kg, 10 kg (4x 2.5 kg), 25 kg
Cat. No. 1.15161	Histosec™ pastilles (without DMSO) solidification point 56-58°C embedding agent for histology	10 kg (4x 2.5 kg), 25 kg

Hazard classification

Cat. No. 1.09204

Please observe the hazard classification printed on the label and the information given in the safety data sheet.

The safety data sheet is available on the website and on request.

Main components of the product

Cat. No. 1.09204		
C.I.52015 + Azure		4.1 g/l
C.I.45380		2.4 g/l
contains CH ₃ OH		
1 l = 0.99 kg		

Other IVD products

Cat. No. 1.00869	Entellan™ new for cover slipper for microscopy	500 ml
Cat. No. 1.01383	Wright's eosin methylene blue solution for microscopy	100 ml, 500 ml, 2.5 l
Cat. No. 1.02439	Eosin Y-solution 0.5%, alcoholic for microscopy	500 ml, 2.5 l
Cat. No. 1.03999	Formaldehyde solution min. 37% free from acid stabilized with about 10% methanol and calcium carbonate for histology	1 l, 2.5 l, 25 l
Cat. No. 1.05174	Hematoxylin solution modified acc. to Gill III for microscopy	500 ml, 1 l, 2.5 l
Cat. No. 1.05175	Hematoxylin solution modified acc. to Gill II for microscopy	500 ml, 2.5 l
Cat. No. 1.05387	Leishman's eosin methylene blue solution modified for microscopy	500 ml
Cat. No. 1.09844	Eosin Y-solution 0.5% aqueous for microscopy	1 l, 2.5 l
Cat. No. 1.11661	Hemacolor® Rapid staining of blood smear staining kit for microscopy	1 set
Cat. No. 1.17081	Eosin Y solution 1%, alcoholic for microscopy	1 l

General remark

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

Literature

- Atlas der klinischen Hämatologie, Löffler, Rastetter, Haferlach, 2004, Springer Verlag 6. Auflage
- Theory and Practice of Histological Techniques, John d. Bancroft, Marilyn Gamble, 2008, Churchill Livingstone Elsevier, Sixth Edition
- Histological & Histochemical Methods, J. A. Kiernan, 1990, Pergamon Press, Second Edition
- Romeis - Mikroskopische Technik, Editors: Maria Mulisch, Ulrich Welsch, 2015, Springer Spektrum, 19. Auflage
- Sobotta, Lehrbuch Histologie, Welsch, 2006, Urban&Fischer, 2. Auflage
- Atlas of clinical cytology, Paul Lopez Cardozo, EM editon medicin
- Klinische Hämatologie, Herbert Begemann, 1975, Georg Thieme Verlag, 2. Auflage
- Conn's Biological Stains: A Handbook of Dyes, Stains and Fluorochromes for Use in Biology and Medicine, 10th Edition, (ed. Horobin, R.W. and Kiernan, J.A). Bios, 2002

H225: Highly flammable liquid and vapor.

H301 + H311 + H331: Toxic if swallowed, in contact with skin or if inhaled.

H317: May cause an allergic skin reaction.

H370: Causes damage to organs (Eyes, Central nervous system).

P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P233: Keep container tightly closed.

P280: Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P310: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P303 + P361 + P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P304 + P340 + P311: IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

Consult instructions
for use

Manufacturer

Catalog number

Batch code

Caution, consult
accompanying documents

Use by
YYYY-MM-DD

Temperature
limitation

Status: 2022-Oct-04

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly available resources.

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