Lyphochek® Quantitative Urine Control Normal (1) and Abnormal (2)

REF 376 Normal (1) 12 x 10 mL
377 Abnormal (2) 12 x 10 mL
375X Bilevel MiniPak 2 x 10 mL

EXP 2023-06-30

63470 Level 1 63471
Level 2 63472

ENGLISH

INTENDED USE

Lyphochek Quantitative Urine Control is intended for use as an assayed quality control urine to monitor the precision of laboratory testing procedures for the analytes listed in this package insert.

SUMMARY AND PRINCIPLE

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices. Two levels of control are available to allow performance monitoring within the clinical range.

For customers in Germany: Quality control materials are required for assessment of laboratory performance as described in the "Guideline for Quality Assurance of Medical Laboratory Examinations following the German Medical Association" (Rili-BÄK regulation).

REAGENT

This product is prepared from human urine with added constituents of human and animal origin, therapeutic drugs, chemicals, hormones and stabilizers. The control is provided in lyophilized form for increased stability.

STORAGE AND STABILITY

This product will be stable until the expiration date when stored unopened at 2 to 8°C. Once the control is reconstituted, all analytes will be stable for 5 days when stored tightly capped at 2 to 8°C, with the following exceptions: Coproporphyrin, Homovanillic Acid (HVA), Dopamine, Epinephrine, Norepinephrine, Metanephrines, Total Porphyrins, Uroporphyrin, 5-Aminolevulinic Acid and Vanillylmandelic Acid (VMA) must be assayed immediately following reconstitution.

The stability of 5-Hydroxyindole Acetic Acid (5-HIAA), Dopamine, Epinephrine, Norepinephrine, Metanephrines, Homovanillic Acid (HVA), 5-Aminolevulinic Acid and Vanillylmandelic Acid (VMA) may be extended to 30 days when the procedures for acidification are followed as described below and the vial is stored tightly capped at 2 to 8°C. After acidification, the control should only be used for the corresponding analytes.

- For 5-Hydroxyindole Acetic Acid (5-HIAA) and 5-Aminolevulinic Acid acidify with 100.0 µL glacial acetic acid per 10.0 mL of reconstituted control.
- For Metanephrines, Dopamine, Epinephrine, Norepinephrine, Homovanillic Acid (HVA) and Vanillylmandelic Acid (VMA), acidify with 80.0 µL 6M HCl per 10.0 mL of reconstituted control or reconstitute the lyophilized control with 10 mL of 0.05M HCl.

This product is shipped under refrigerated conditions.

RECONSTITUTION

Using a volumetric pipet, reconstitute each vial with 10.0 mL of distilled or deionized water. Replace the stopper and allow the control to stand for at least 15 minutes, swirling occasionally. Before sampling, gently invert the vial several times to ensure homogeneity.

PROCEDURE

This product should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit, or reagent being used.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities. In the event of damage to packaging, contact the local Bio-Rad Laboratories Sales Office or Bio-Rad Laboratories Technical Services.

LIMITATIONS

- This product should not be used past the expiration date.
- If there is evidence of microbial contamination or excessive turbidity in the reconstituted control, discard the vial.
- This product is not intended for use as a standard.

ASSIGNMENT OF VALUES

The mean values and the corresponding $\pm 3SD$ ranges printed in this insert were derived from replicate analyses and are specific for this lot of product. Data from Unity™ Interlaboratory Program are included in the determination of some ranges. The tests listed were performed by the manufacturer and/or independent laboratories using manufacturer supported reagents and a representative sampling of this lot of product. It is recommended that each laboratory establish its own acceptable ranges and use those provided only as guides. Laboratory established ranges may vary from those listed during the life of this control. [Customers in Germany have to follow the requirements as described in the Rili-BÄK regulation.] Variations over time and between laboratories may be caused by differences in laboratory technique, instrumentation and reagents, or by manufacturer test method modifications.

Unity™ Interlaboratory Program is a proprietary Bio-Rad software program holding more than 2 billion QC data points from thousands of laboratories.

Refer to www.qcnet.com for insert update information.

SPECIFIC PERFORMANCE CHARACTERISTICS

This product is a freeze-dried product manufactured under rigid quality control standards. To obtain consistent vial-to-vial assay values, the control requires proper storage and handling as described.

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ENGLISH

Biological source material. Treat as potentially infectious.

The serum from each donor contributing urine for this product was tested as required by FDA accepted methods. Tests results were non-reactive or negative for evidence of infection due to Human Immunodeficiency Virus (HIV), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV). This product may also contain other human source materials for which there are no approved tests. In accordance with good laboratory practice, all human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

[Level Abnormal](#)

Warning

Hazard (H) and Precautionary (P) Statements

H315 Causes skin irritation. **H319** Causes serious eye irritation. **P280** Wear protective gloves / protective clothing / eye protection / face protection. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P363** Wash contaminated clothing before reuse. **P333+P313** If skin irritation or rash occurs: Get medical advice / attention. **P337+P313** If eye irritation persists: Get medical advice / attention. **P302+P352** IF ON SKIN: Wash with plenty of soap and water.

Safety Data Sheet (SDS) available for professional users on [www.bio-rad.com](#).

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