Third Generation TSH

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE[®] 2000

Third Generation TSH

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE[®] 2000 Systems Analyzers — for the quantitative measurement of thyrotropin (TSH) in serum. Measurements of thyroid stimulating hormone produced by the anterior pituitary are used in the diagnosis of thyroid or pituitary disorders.

Catalog Numbers:
L2KTS2 (200 tests)
L2KTS6 (600 tests)

Test Code: **TSH** Color: **Red**

Summary and Explanation

Thyroid stimulating hormone (thyrotropin, TSH) is a pituitary hormone which, through its action on the thyroid gland, plays a major role in maintaining normal circulating levels of the iodothyronines, T4 and T3. TSH is controlled by negative feedback from circulating T4 and T3, and by the hypothalamic hormone TRH (thyrotropin releasing hormone). TSH exhibits a small circadian rhythm.

In primary hypothyroidism, where there is impaired production of thyroid hormones, the TSH level is typically highly elevated. In secondary or tertiary hypothyroidism, on the other hand, where thyroid hormone production is low as a consequence of pituitary or hypothalamic lesions, the TSH level is usually low. In hyperthyroidism, the TSH level is typically suppressed to subnormal levels. Less often, this condition may result from hyperstimulation of the thyroid, due to hypothalamic or pituitary lesions, in which case the TSH level is usually increased.

Measurement of circulating TSH has been used as a primary test for differential diagnosis of hypothyroidism and as an aid in monitoring the adequacy of thyroid hormone replacement therapy. It should be remembered that hyperthyroidism and hypothyroidism are graded conditions. This implies that not all patients in these disease categories can be expected to

have TSH levels far outside the euthyroid range. On the other hand, TSH levels exit the euthyroid reference range in the very early phases of developing thyroid disease, while the patient's disease is still subclinical and thyroid hormone levels remain within their euthyroid reference ranges.

Research studies have found that the apparently healthy patients with TSH > 2.0 $\mu\text{IU/mL}$ have increased risk to develop thyroid diseases in the next 20 years. It has been suggested that it is likely that the upper limit of the serum TSH euthyroid reference range will be reduced to 2.5 $\mu\text{IU/mL}$ because > 95% of rigorously screened normal euthyroid volunteers have serum TSH values between 0.4 and 2.5 $\mu\text{IU/mL}$.⁷

Principle of the Procedure

IMMULITE 2000 Third Generation TSH is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 60 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

EDTA plasma should not be used as a sample type.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot

activators and/or anticoagulants. IMMULITE 2000 Third Generation TSH has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 75 µL serum

Storage: 5 days at 2–8°C, or 1 month at –20°C.¹

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.⁸⁻¹⁰

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412 P273, P501	Harmful to aquatic life with long lasting effects. Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations. Contains: 2-methyl-2H-isothiazol-3-one; TSH Adjustors
----------------------------------	--

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

TSH Bead Pack (L2TS12)

With barcode. 200 beads, coated with monoclonal murine anti-TSH antibody. Stable at 2–8°C until expiration date.

L2KTS2: 1 pack
L2KTS6: 3 packs

TSH Reagent Wedge (L2TSA2)

With barcode. 23 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal goat anti-TSH antibody in buffer (with preservative), dispensed equally into chambers A and B. Stable at 2–8°C until expiration date.

L2KTS2: 1 wedge
L2KTS6: 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

TSH Adjustors (LTSL, LTSH)

Two vials (Low and High) of lyophilized human TSH in a serum/buffer matrix. Reconstitute each vial with **4.0 mL** distilled or deionized water. Let stand for 30 minutes, then mix by *gentle* swirling or inversion. Stable at 2–8°C for 30 days, after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KTS2: 1 set
L2KTS6: 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

TSH Sample Diluent (L2TSZ)

For on-board dilution of high samples. 25 mL of concentrated (ready-to-use) TSH-free serum/buffer matrix. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2TSZ: 3 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LTGCM: Third Generation TSH Control (single-level, low control)

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:

4 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of TSH.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Expected Values

Based on its relationship to IMMULITE Third Generation TSH (see Method Comparison), the assay can be expected to have essentially the same euthyroid reference range:

Adult Euthyroid: 0.4–4 µIU/mL

TSH levels in frank hyperthyroidism are typically less than 0.01 µIU/mL. Values between 0.01 µIU/mL and 0.4 µIU/mL should be further evaluated, as they may indicate borderline hyperthyroidism, or may be due to severe nonthyroidal illness (NTI) or acute drug therapy.¹

Pediatric: Reference intervals for the pediatric population (children and adolescents) were established for the IMMULITE Third Generation TSH assay in accordance with CLSI guideline EP28-A3C.¹¹

For analysis of data, the population was divided into three age subgroups:

- Infants: subjects aged 1–23 months
- Children: subjects aged 2–12 years
- Adolescents: subjects aged 13–20 years

A non-parametric approach was used to establish the reference intervals for children and adolescents where the 2.5 and 97.5 percentiles of the distribution of values were calculated. For the infant population, a robust measure of location and spread, as developed by Horn and Pesce, was used to estimate the 2.5 and the 97.5 percentile reference intervals, accommodating the smaller sample size.^{11–13}

The reference intervals detailed by age group and number of samples are presented in the Reference Intervals table.

IMMULITE 2000 Third Generation TSH Pediatric Reference Intervals

Age Group	n	Conventional (μ IU/mL)	SI (mIU/L)
Infants* (1 – 23 Months)	90	0.83 – 6.5	0.83 – 6.5
Children (2 – 12 Years)	195	0.58 – 4.1	0.58 – 4.1
Adolescents (13 – 20 Years)	148	0.39 – 4.0	0.39 – 4.0

* The upper limit (97.5 percentile) of the infant reference interval was determined to be 6.5 μ IU/mL (mIU/L). Data from this infant population (n = 90) have demonstrated a highly skewed distribution to the right; therefore, the estimate of the upper limit of the reference interval has some uncertainty with a 90% probability that the upper limit of the reference interval can be between 5.6 to 7.7 μ IU/mL (mIU/L).

Consider these limits as guidelines only. Each laboratory should establish its own reference ranges.

Limitations

The IMMULITE 2000 Third Generation TSH assay is not indicated for use with neonatal blood spots for newborn screening.

As with any immuno-recognition measurement of a peptide, extremely rare genetic variants may exhibit varying degrees of detection.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data representative of the assay's performance. Results are expressed in μ IU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Calibration Range: Up to 75 μ IU/mL (WHO 2nd IRP 80/558)

Analytical Sensitivity: 0.004 μ IU/mL

High-Dose Hook Effect: None up to 14,000 μ IU/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three TSH solutions (273, 1387 and 3290 μ IU/mL) were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for TSH. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of packed red blood cells in concentrations up to 30 μ L/mL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 2000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 10 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. All samples were assayed by the IMMULITE 2000 Third Generation TSH procedure. By linear regression:

(Heparin) = 0.95 (Serum) + 0.06 μ IU/mL
r = 0.999

(EDTA) = 0.64 (Serum) + 0.02 μ IU/mL
r = 0.998

(SST) = 0.99 (Plain Tubes) – 0.11 μ IU/mL
r = 0.991

Means:

1.65 μ IU/mL (Serum)
1.62 μ IU/mL (Heparin)
1.07 μ IU/mL (EDTA)
1.52 μ IU/mL (SST)

Method Comparison: The assay was compared to IMMULITE Third Generation TSH on 106 samples. (Concentration range: approximately 0.006 to 60 μ IU/mL. See graph.) By linear regression:

(IML 2000) = 0.94 (IML) + 0.23 μ IU/mL
r = 0.985

Means:

6.4 μ IU/mL (IMMULITE 2000)
6.5 μ IU/mL (IMMULITE)

References

- 1) Hay ID, Bayer MF, et al. American Thyroid Association assessment of current free thyroid hormone and thyrotropin measurements and guidelines for future clinical assays. *Clin Chem* 1991;37:2002-8.
- 2) Lindstedt G, et al. Clinical use of laboratory thyroid tests and investigations. *JIFCC* 1994;6:136-41.
- 3) Mandel SJ, Brent GA, Larsen PR. Levothyroxine therapy in patients with thyroid disease. *Ann Intern Med* 1993;119:492-502.
- 4) Singer PA, Cooper DS, et al. Treatment guidelines for patients with hyperthyroidism and hypothyroidism. *JAMA* 1995;273:808-12.
- 5) Spencer CA, Takeuchi M, Kazarosyan M, et al. Interlaboratory/intermethod differences in functional sensitivity of immunometric assays of thyrotropin (TSH) and impact on reliability of measurement of subnormal concentrations of TSH. *Clin Chem* 1995;41:367-74.
- 6) Spencer CA, Takeuchi M, Kazarosyan M. Current status and performance goals for serum thyrotropin (TSH) assays. *Clin Chem* 1996;42:140-5.
- 7) Baloch Z, Carayon P, Conte-Devolx B, Demers LM, Feldt-Rasmussen U, Henry JF, et al.; Guidelines Committee, National Academy of Clinical Biochemistry (NACB). Laboratory medicine practice guidelines (LMPG). Laboratory support for the diagnosis and monitoring of thyroid disease. *Thyroid* 2003 Jan;13(1):3-126. Also available at http://www.nacb.org/lmpg/thyroid_LMPG_Word.stm (accessed January 2005).
- 8) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. *MMWR*, 1988;37:377–82, 387–8.
- 9) Clinical and Laboratory Standards Institute

(formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 10) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030. 11) Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Guideline EP28-A3C. 12) Horn PS, Pesce AJ, Reference Intervals. A User's Guide, Washington, DC:AACC Press; 2005. 13) Reed AH, Henry RJ, Mason WB. Influence of statistical method used on the resulting estimate of normal range. 1971; 17:275-284.

Technical Assistance

Contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (μ IU/mL)

	Mean ³	Within-Run ¹		Total ²	
		SD ⁴	CV ⁵	SD	CV
1	0.016	0.002	12.5%	0.002	12.5%
2	0.32	0.017	5.3%	0.017	5.3%
3	1.3	0.05	3.8%	0.06	4.6%
4	3.3	0.13	3.9%	0.16	4.8%
5	7.3	0.37	5.1%	0.37	5.1%
6	19	0.73	3.8%	0.85	4.5%
7	39	2.0	5.1%	2.5	6.4%

Linearity (µIU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	16 in 16 ⁵	6.4	—	—
	8 in 16	3.1	3.2	97%
	4 in 16	1.6	1.6	100%
	2 in 16	0.86	0.80	108%
	1 in 16	0.42	0.40	105%
2	16 in 16	21	—	—
	8 in 16	10.3	10.6	97%
	4 in 16	4.7	5.3	89%
	2 in 16	2.4	2.6	92%
	1 in 16	1.1	1.3	85%
3	16 in 16	58	—	—
	8 in 16	30	29	103%
	4 in 16	15	15	100%
	2 in 16	7.6	7.3	104%
	1 in 16	3.7	3.6	103%

Recovery (µIU/mL)


	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	0.78	—	—
	A	14	14	100%
	B	71	70	101%
	C	166	165	101%
2	—	4.7	—	—
	A	18	18	100%
	B	75	74	101%
	C	164	169	97%
3	—	7.2	—	—
	A	20	20	100%
	B	78	76	103%
	C	173	171	101%

Specificity

Compound ¹	ng/mL Added ²	% Cross reactivity ³
FSH	100	ND
LH	200	ND
HCG	1000	ND

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2018-2023 Siemens Healthcare Diagnostics. All rights reserved.
ADJUSTOR AB