



# Specification

---

1.04699.0100 Immersion oil for microscopy

---

---

## Specification

---

Refractive index (n <sub>20/D</sub> )	1.515 - 1.517	
Density (d <sub>20 °C/4 °C</sub> )	1.0245 - 1.0265	
Transmission (380 nm; 1 cm)	≥ 65	%
Transmission (400 nm; 1 cm)	≥ 78	%
Transmission (450 nm; 1 cm)	≥ 90	%
Fluorescence (as quinine at 365 nm)	≤ 1500	ppb
Viscosity (20 °C)	100 - 120	mPa*s
Suitability for microscopy	passes test	

Dr. Karl-August Reiffen

---

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.